

The diagram illustrates the relationship between Microscope Specifications and Image Acquisition Settings. On the left, a box labeled "MICROSCOPE Specifications" is connected to a "Children Entities (choose one)" container. This container lists "Arc" and "Filament" under the "Type" category. Arrows point from these entities to the "Image Acquisition Settings" box on the right. The "Image Acquisition Settings" box contains two "WavelengthRange" sections. The first section, associated with "Arc", lists "PeakWavelength", "CutOn", "CutOff", and "WavelengthProfileFile". The second section, associated with "Filament", lists "WavelengthRange".

